


OFFLOAD CANADA INC.

56 Aberfoyle Crescent
Toronto ON M8X 2W4
www.offloadcanada.com
info@offloadcanada.com
+1 (416) 729-6943



ITEM.1.4



ITEM.1.4
WALL SUPPORT FOR HANDS GEL.
MADE ENTIRELY OF STAINLESS STEEL

OFFLOAD CANADA INC.

56 Aberfoyle Crescent
Toronto ON M8X 2W4
www.offloadcanada.com
info@offloadcanada.com
+1 (416) 729-6943



ITEM.1.5

	<p>Automatic measuring devices with 316 stainless steel temperature gauge precision +/- 0.2. Power supply batteries stylus or powerbank included. Detection speed 1 sec It can work without the help of a supervisor.</p> <p>This detector has proven extremely effective in large inflows, works either with the wrist or with the forehead and à a tolerance of +o- 0.2 degrees. In case the person is not able to carry out the measurement with the wrist, it is enough to bring the forehead about 10 cm from the device and the detection will be performed with the same precision detection time 1 second. The stand that supports it has been designed with an inclination that allows the taller people to measure with the forehead, and at the same time ensure even the lower people the possibility of doing so. An example of our installation are schools where with 3 stands where 3400 students can measure in 20 minutes. The thermo scanner should be placed inside and not below 8 degrees.</p>
Weight	10 kg
Height	120cm
Width	20 cm
Depth	30 cm
Power supply	works with 4 AA batteries or powerbank DC 5V 500Ma
Temperature detection speed	1 secondo
Audible alarm	37.5 Audible and visual alarm
Operating temperature	10 'C / 35 'C
Outputs controlled	Yes

Distributed by Offload Canada Inc.

Te No. +1 (416) 729-6943

www.offloadcanada.com